



---

# News Release

---

FOR IMMEDIATE RELEASE

30 October, 2018

LARRY VAN BOGAERT

P.O. BOX 585  
Elkhorn, NE 68022-0585

## NEW PRODUCT RELEASE

SP200VCA Variable Flow Pump

877.324.5444 Toll-Free  
217.498.0342 TEL  
217.498.0352 FAX

### **Low-cost, Variable-flow Peristaltic Pump**

Compact pump provides variable-flow and exceptional life performance in a small attractive package. Enjoy the self-priming, run-dry, reversible, fluid-only-in-tubing benefits of a peristaltic pump in your fluid transfer or metering application. Flow rates are available up to 320 ml/min (5.07 gph) with a choice of motor speeds and tubing sizes. Select from a wide range of tubing materials for sterile operation, general purpose transfer or chemical handling. Splash-resistant housing mounts easily on the wall or on the bench. Manual and remote control options.

A long time producer and supplier of fluid handling pumps and systems, APT Instruments offers a wide range of small, compact solutions for fluid and air handling. In addition, application engineering assistance is provided to meet special product or custom design requirements. OEM and Distributor Programs. – APT Instruments, P.O. Box 585, Elkhorn, NE 68022-0585, 877.324.5444, 217.498.0342, FAX 217.498.0352.

[sales@aptinstruments.com](mailto:sales@aptinstruments.com), [www.aptinstruments.com](http://www.aptinstruments.com).

END

Reference photo included – contact APT for HiRes.